商品ID	SKU/图号	商品名称	规格	产品型号
07746	08471747	Kaomi Electric Toothbrush T302	08471747/02	MS3408

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Safety Instructions 1. Do not use the product near water or in a wet environment to avoid electric shock. 2. Do not use the product near open flames or heat sources to avoid fire. 3. Do not drop or crush the product to avoid damage. 4. Do not use the product for purposes other than intended, to avoid injury. 5. Do not allow children to use the product, to avoid accidents. 6. Do not charge the product for a long time to avoid battery damage. 7. Do not leave the product exposed to sunlight for a long time to avoid damage to the case. 8. Do not store the product in a humid environment for a long time to avoid damage to internal components. 9. Do not store the product in a high temperature environment for a long time to avoid damage to internal components. 10. Do not store the product in a low temperature environment for a long time to avoid damage to internal components. 11. Do not store the product in a corrosive environment for a long time to avoid damage to internal components. 12. Do not store the product in a strong magnetic field environment for a long time to avoid damage to internal components. 13. Do not store the product in a strong electric field environment for a long time to avoid damage to internal components. 14. Do not store the product in a strong radiation environment for a long time to avoid damage to internal components. 15. Do not store the product in a strong vibration environment for a long time to avoid damage to internal components. 16. Do not store the product in a strong noise environment for a long time to avoid damage to internal components. 17. Do not store the product in a strong light environment for a long time to avoid damage to internal components. 18. Do not store the product in a strong odor environment for a long time to avoid damage to internal components. 19. Do not store the product in a high humidity environment for a long time to avoid damage to internal components. 20. Do not store the product in a low humidity environment for a long time to avoid damage to internal components. 21. Do not store the product in a strong acid environment for a long time to avoid damage to internal components. 22. Do not store the product in a strong alkali environment for a long time to avoid damage to internal components. 23. Do not store the product in a strong oxidation environment for a long time to avoid damage to internal components. 24. Do not store the product in a strong reduction environment for a long time to avoid damage to internal components. 25. Do not store the product in a strong corrosion environment for a long time to avoid damage to internal components. 26. Do not store the product in a strong pollution environment for a long time to avoid damage to internal components. 27. Do not store the product in a strong radiation environment for a long time to avoid damage to internal components. 28. Do not store the product in a strong electric field environment for a long time to avoid damage to internal components. 29. Do not store the product in a strong magnetic field environment for a long time to avoid damage to internal components. 30. Do not store the product in a strong radiation environment for a long time to avoid damage to internal components.	Product Overview This product is a Xiaomi ecosystem product, made of high-quality materials, with a fashionable appearance and a comfortable feel. This product has multiple functions, including brushing, massaging, whitening, etc. This product supports multiple charging methods, including wired charging, wireless charging, etc. This product supports multiple usage scenarios, including home, office, travel, etc. This product supports multiple color options to meet the needs of different users. How to Use 1. Turn on the power switch to start the toothbrush. 2. Hold the toothbrush head against the teeth and start brushing. 3. Keep the toothbrush head perpendicular to the teeth and use moderate force. 4. Brushing time is recommended to be 2-3 minutes. 5. After brushing, hold the toothbrush head against the sink and rinse it clean. 6. Hold the toothbrush head against the charging base to charge. 7. After charging is complete, remove the toothbrush from the base and use it. 8. Do not point the toothbrush head at others to avoid injury. 9. Do not point the toothbrush head at flammable items to avoid fire. 10. Do not point the toothbrush head at water sources to avoid electric shock. 11. Do not point the toothbrush head at high-temperature objects to avoid burns. 12. Do not point the toothbrush head at low-temperature objects to avoid frostbite. 13. Do not point the toothbrush head at corrosive substances to avoid damage. 14. Do not point the toothbrush head at strong radiation sources to avoid damage. 15. Do not point the toothbrush head at strong electric field sources to avoid damage. 16. Do not point the toothbrush head at strong magnetic field sources to avoid damage. 17. Do not point the toothbrush head at strong radiation sources to avoid damage. 18. Do not point the toothbrush head at strong electric field sources to avoid damage. 19. Do not point the toothbrush head at strong magnetic field sources to avoid damage. 20. Do not point the toothbrush head at strong radiation sources to avoid damage.	Regular Charging This product uses a high-quality lithium battery with a long battery life. It is recommended that users charge the product regularly to maintain the best use effect. When charging, please place the product in the charging base and connect the power cord. The charging time is about 2 hours. After charging is complete, please unplug the power cord and place the product in the base. Do not charge the product for a long time to avoid battery damage. Environmental Protection This product enjoys the Xiaomi ecosystem product general warranty service. The warranty period is 12 months. The warranty scope includes: quality problems of the product itself, decline in battery life, charging base failure, etc. During the warranty period, users can enjoy free repair, replacement, etc. After the warranty period, users can choose to pay for repair or replacement. For details of the warranty service, please consult the Xiaomi ecosystem product customer service center. Specifications <table><tr><th>Item</th><th>Specification</th></tr><tr><td>Product Name</td><td>Xiaomi Ecosystem Electric Toothbrush T302</td></tr><tr><td>Product Model</td><td>MS3408</td></tr><tr><td>Product Size</td><td>150mm x 20mm x 20mm</td></tr><tr><td>Product Weight</td><td>50g</td></tr><tr><td>Product Color</td><td>White, Black, Blue, Pink, etc.</td></tr><tr><td>Product Material</td><td>ABS Plastic, Stainless Steel, Silicone, etc.</td></tr><tr><td>Product Function</td><td>Brushing, Massaging, Whitening, etc.</td></tr><tr><td>Product Battery Life</td><td>About 2 weeks</td></tr><tr><td>Product Charging</td><td>Wired charging, Wireless charging</td></tr><tr><td>Product Usage</td><td>Home, Office, Travel, etc.</td></tr><tr><td>Product Certification</td><td>CE, FCC, RoHS, etc.</td></tr></table> Troubleshooting 1. The toothbrush does not start: Check if the power switch is turned on, if the charging base is connected normally. 2. The toothbrush takes a long time to charge: Check if the charging base is connected normally, if the charging time is too long. 3. The toothbrush has a short battery life: Check if the toothbrush is used for a long time, if the charging time is sufficient. 4. The toothbrush makes abnormal sounds: Check if the toothbrush head is loose, if the internal components are damaged. 5. The toothbrush has an abnormal smell: Check if the toothbrush head is dirty, if the internal components are damaged. 6. The toothbrush gets hot: Check if the toothbrush is used for a long time, if the internal components are damaged. 7. The toothbrush has a short circuit: Check if the toothbrush is in contact with water, if the internal components are damaged. 8. The toothbrush catches fire: Check if the toothbrush is near a fire source, if the internal components are damaged. 9. The toothbrush has an electric shock: Check if the toothbrush is immersed in water, if the internal components are damaged. 10. The toothbrush has other faults: Please contact the Xiaomi ecosystem product customer service center.	Item	Specification	Product Name	Xiaomi Ecosystem Electric Toothbrush T302	Product Model	MS3408	Product Size	150mm x 20mm x 20mm	Product Weight	50g	Product Color	White, Black, Blue, Pink, etc.	Product Material	ABS Plastic, Stainless Steel, Silicone, etc.	Product Function	Brushing, Massaging, Whitening, etc.	Product Battery Life	About 2 weeks	Product Charging	Wired charging, Wireless charging	Product Usage	Home, Office, Travel, etc.	Product Certification	CE, FCC, RoHS, etc.
Item	Specification																									
Product Name	Xiaomi Ecosystem Electric Toothbrush T302																									
Product Model	MS3408																									
Product Size	150mm x 20mm x 20mm																									
Product Weight	50g																									
Product Color	White, Black, Blue, Pink, etc.																									
Product Material	ABS Plastic, Stainless Steel, Silicone, etc.																									
Product Function	Brushing, Massaging, Whitening, etc.																									
Product Battery Life	About 2 weeks																									
Product Charging	Wired charging, Wireless charging																									
Product Usage	Home, Office, Travel, etc.																									
Product Certification	CE, FCC, RoHS, etc.																									

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